

Cleaning Mask

1. Hand wash mask in warm water with a mild dishwashing detergent.
DO NOT use bleach, alcohol or cleaning solutions containing alcohol.
DO NOT use cleaners containing conditioners or moisturizers.
2. Rinse thoroughly. Pat dry with a towel and air dry completely before use.
3. Inspect all parts for damage or wear: replace if necessary.

Cleaning the Headgear

The headgear should be cleaned weekly by hand or machine wash in the gently cycle. Lay the headgear flat or line dry. **DO NOT** place the headgear in the dryer.